



Product Name : " Blood Collection Monitor "  
Product Code : " SI 31 "



[chhabrainstruments.com](http://chhabrainstruments.com)

Description :

### **Blood Collection Monitor:-**

We are prominent suppliers of Blood Collection Monitors, which are compact instruments, providing smooth & gentle rocking for homogeneous mixing of blood with anticoagulant. The entire process happens without clotting of blood cells during collection of blood from donor. These fully automatic microprocessor based monitor machines have the capacity of mixing blood at a rate of 10-12 cycles/min. When the blood volume reaches preset value, the clamp starts working and the process is automatically terminated. The LCD Screen of the monitor shows the volume of collected blood and time taken in the process. These machines have the battery back -up of 4-20 hours, depending upon the requirement of the task.