

Product Name : " Blood Shaker / Mixer "
Product Code : " SI 34 "



Description :

Blood Shaker / Mixer:-

We are manufacturers and suppliers of an extensive range of Blood Mixers, which is used by Hospitals, Pathological, Bio-medical Research and Institutional Laboratories. Blood Mixers are compactly designed equipment to provide smooth, intermittent and gentle rocking for homogeneous mixing of blood with anticoagulant. These machines are Instrumental in prevention of blood clotting during collection of blood from donor. Also, these machines are easy to operate and maintain.