

Product Name : " Donor Couch "
Product Code : " SI 37 "



chhabrainstruments.com

Description :

Donor Couch:-

We are engaged in offering **Donor Couch**, which is a completely upholstered donor chair designed to provide comfort and safety to the donors.

- Electric Donor Couch is designed to provide a comfort postman for the Donor.
- Donor's position can be easily shifted by micro motion at your whole at the touch of a button and remote it is easy to operate for all with three types of position.
- All surface are clean & smooth, easily accessible and can be kept clean.
- Electrical & mechanical accessories are housed inside the structure making the product safer.
- Tilt adjustments can be done by using remote control, ensuring safety & comfort to the donor.
- Caster wheels with locking system.
- Inbuilt heavy duty geared motor.
- We make available uniquely designed Blood Donor Couch at market leading prices.

Technical Specification:

Power Supply	220-240 VAC
Operating Voltage	24 VDC
Weight	150 kg approx.
Lifting Capacity	up to 150 kg
Control	Remote
Overall Length	1900 mm
Seat Width	550 mm
Overall Width	660 mm
Seat Height	550 mm