

Product Name : " VDRL Shaker "

Product Code : " SI 39 "



Description :

VDRL Shaker:-

- Compact shaker for mixing blood samples.
- Timer options for accurate mixing time.
- Digital RPM controller and digital timer available for exact results.
- Stable at high speeds.
- Wide speed range.
- Separate platform attachment available for mixing more than one samples.

Technical Data

Type of movement	Rotational
Speed range	0- 350 RPM
Operational mode	Continuous operation
Timer option	Yes
Dimensions	37 x 37 x 18

Weight	11 kgs
Permissible ambient temperature	5 - 40 deg.
Permissible relative moisture	80%
Motor type	PMDC motor
Maximum load	6- kg
Speed display	None
Supply	220 V, 50 Hz.