

Product Name : " TOOL SET "
Product Code : " R.A.C 27 "



Description :

TOOL SET:-

AIM:-

Tools for the assembly of the frame work, the hydraulics (piping) and electrics (wiring)

- Leak detection in a refrigeration system

TECHNICAL SPECIFICATION:-

Consists of a filled tool box made of plastic. The following tools will be supplied:-

- 1 No. : Pipe cutter
- 1 No. : Adjustable spanner
- 1 No. : Crimping tool with threaded cutter for bolts,

-
- 1 No. : Tube flaring unit
- 1 No. : Side cutter
- 1 No. : Cable stripping knife
- 1 No. : Crimping tool for wire ends
- 1 No. : Wire stripper
- 1 Set. : Set of pipe bending pliers (3 pieces)
- 1 Set. : Set of screw drivers, completely insulated (4x slotted, 2x Philips)
- 1 Set. : Set of combination spanners (19 pieces)
- 1 No. : A leakage detector for several refrigerant gases and a voltage tester.

Specification

1. Commercially available tools for assembly of refrigeration systems.
2. Battery-operated leak detector suitable for refrigerant R134a.
3. Voltage tester for direct current and alternating current
Voltage tester: 12...400VAC/DC