

Product Name : " RECIPROCATING PUMP TEST SETUP "
Product Code : " R.A.C 47 "



Description :

RECIPROCATING PUMP TEST SETUP:-

DESCRIPTION:-

[Best Cheap UK Replica Watches For Sale](#)

The Set-up is designed to study the performance of Reciprocating Pump. The Set-up consists of a Reciprocating Pump coupled with electrical motor, supply tank, measuring tank & pipe fittings for closed loop water circulation. Pressure and Vacuum Gauge are connected on delivery and suction side of pump for the purpose of measurement. The flow rate of water is measured using measuring tank and stop watch provided.

FEATURES:-

- Closed loop water circulation
- Compact & stand alone set-up
- MS Excel sample calculation program on demand
- Stainless Steel tanks and wetted parts
- Superb painted structure

- Simple to Operate & Maintain

TECHNICALS SPECIFICATIONS	DESCRIPTION
Pump Best UK Swiss Replica Watches	Double acting, single Cylinder, Capacity 1 Head 5 Kg/cm ² (max.)
Drive	AC Motor with step cone pulley arrangement Motor with DC Drive along with non-contact
Supply Tank	Capacity 70 Ltrs
Measuring Tank	Capacity 50 Ltrs. fitted with Piezometer Tu
Piping	GI / PVC
Stop Watch	Electronic
Control Panel	With required electrical instrumentation