



Product Name : " WATER COOLER "
Product Code : " CKE 32 "



Description :

WATER COOLER:-

Outer body & inner chamber of stainless steel sheet. We are the prominent Manufacture and Supplier of premium quality Drinking Water Cooler which we supply all over India to cater the needs of our clients. Our Drinking Water Cooler is available at the affordable market price for our customers. The product is specifically made to provide pure, safe and cool drinking water. Water Cooler Inner Chamber 304grade & Outer Body 202 grade of Stainless Steel Sheet. Puff Insulation, 2 Outlet tab of water.

Water Cooler Made by	Chhabra Engineers & Instruments
Compressor	Emerson Copland / Tecumseh India
Condenser	Fins and copper Tube Air Cooled

Motor capacity	10/16/25 W
Tank Capacityltrs
Cooling Capacityltrs
Inner Chamber	304 grade
Outer body	202 grade
Refrigerant	
Insulation	Puff Insulation
Cooling pipe	Copper tube (Met tube, Total Line Carrier)
Outlet Tab	25 to 60 Ltrs Single Tab , 80 to 200 Ltrs 2 /3 tab
Water drain	GI Pipe Fitting & Ball valve (brass)