



Product Name : " DEEP FREEZER "  
 Product Code : " CKE 33 "



Description :

## DEEP FREEZER:-

Outer Body Stainless Steel. Inner Chamber of SS Sheet & Hard top.

<b>Deep Freezer Made by</b>	Chhabra Engineers & Instruments
<b>Compressor</b>	Emerson Copland / Tecumseh India
<b>Condenser</b>	Fins and copper Tube Air Cooled
<b>Motor capacity</b>	16/25 W
<b>Tank Capacity</b>	As per customer requirement.
<b>Inner Chamber</b>	304 grade
<b>Outer body</b>	202 grade
<b>Refrigerant</b>	404 A
<b>Insulation</b>	Puff Insulation
<b>Cooling pipe</b>	Copper tube (Met tube, Total Line Carrier )
<b>Water drain</b>	GI Pipe Fitting & Ball valve (brass)