

Product Name : " Laboratory Stirrer "
Product Code : " SI 09 "



Description :

Laboratory Stirrer:-

The motor stirrer having maximum speed up to 4000 r.p.m. Unit is clamped to C.P. rod fitted on a Heavy cast iron base. The height and angle of the stirrer can be adjusted according to the requirement. Variable speed motor with 1/20 H.P. Regulator is fitted with S.S. stirrer rod to give smooth and noiseless operation.

OPTIONAL:-
i). Automatic 0-60 minutes Timer