

**Product Name :** " Oil Bath with Stirrer "

**Product Code :** " SI 12 "



**Description :**

### **Oil Bath with Stirrer:-**

Double Wall Construction. Inner Chamber made of highly polished stainless steel sheet and the outer is made of mild steel duly finished in white stoving enamel/powder coated paint with mat finished colour combinations. Temperature range 50° C to 250° C is controlled by a thermostat. Supplied complete with stirrer, diffuser shelves, cover etc.

<b>Chamber Size in mm</b>	<b>Chamber size in inches Approx.</b>
a). 250x250x250mm	(10"x10"x10")
b). 300x300x250mm	(12"x12"x10")
c). 350x350x250mm	(14"x14"x10")

#### **OPTIONAL:-**

- i). Mercury contact thermometer and relay with accuracy better than 1°C.
- ii). Digital temp. Indicator-cum-controller.

iii). Microprocessor PID digital temp. Indicator-cum-Controller.

iv). Digital temp. Indicator.