

**Product Name :** " Water Bath Incubator Shaker (Metabolic Shaking Incubator) "

**Product Code :** " SI 13 "



**Description :**

### **Water Bath Incubator Shaker (Metabolic Shaking Incubator):-**

Temperature range from room temperature to 95°C. Thermostatic control with an accuracy of  $\pm 1^\circ\text{C}$ . Double walled, inside made of stainless steel and outside mild steel sheet, painted in epoxy powder coating, Oscillating tray is riding on ball rollers. Shaking speed range 40 to 140 r.p.m. Shaking tray can hold either test tubes or flask of 25ml, 50ml, 100ml. To work on 220/230 volts A.C. without R.P.M. Indicator.

<b>Shaking Tray Size (S.S)</b>	<b>Capacity (Approx.)</b>
a) 275 x 275 x 150mm	12 ltrs.
b) 405 x 300 x 150mm	18 ltrs.

#### **OPTIONAL:-**

- |  |
|--|
| i). Digital temp. Indicator.   |
| ii). Digital temp. Indicator –cum-Controller.  |
| iii). Microprocessor PID Digital Temp. Indicator-cum-Controller.                               |
| iv). Shaking tray for 25ml, 50ml, 100ml, 250ml, 500ml, 25mm dia. Test tubes or conical flasks. |

