

Product Name : " Refrigerated Liquid Bath "

Product Code : " S115 "



Description :

Refrigerated Liquid Bath :-

Standard Construction, inner Chamber made out of stainless steel sheet with digital temp. Controller- cum- indicator. Temperature can be controlled from range minus(-)20°C to plus(+)-20°C. & (-)40°C to plus(+)-20°C.

Temp.	Chamber Size:	Capacity (Approx.)
(-)20°C to plus(+)-20°C	450x300x300mm	42 liters
(-)40°C to plus(+)-20°C	450x300x300 mm	42 liters