

**Product Name :** " Refrigerated Liquid Bath "

**Product Code :** " S115 "



**Description :**

### **Refrigerated Liquid Bath :-**

Standard Construction, inner Chamber made out of stainless steel sheet with digital temp. Controller- cum- indicator. Temperature can be controlled from range minus(-)20°C to plus(+)-20°C.& (-)40°C to plus(+)-20°C.

Temp.	Chamber Size:	Capacity (Approx.)
(-)20°C to plus(+)-20°C	450x300x300mm	42 liters
(-)40°C to plus(+)-20°C	450x300x300 mm	42 liters